# BÜCHI 535

Schmelzpunkt  
Melting Point  
Point de fusion

95340, 10 1993, A4 500, FAE

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## 1. Technical data

|  |   |
|--|---|
| Mains voltage  | 110, 220, 240V; 50/60 Hz  |
| Power input  | 200 Watt  |
| Permissible voltage fluctuation                        | +/- 10%   |
| Voltage for heating                                    | mains voltage   |
| Voltage for illumination                               | 12 Volts a.c.   |
| Fan  | 12 Volts d.c.   |
| Stirrer motor  | mains voltage   |
| Oil volume for melting and boiling point determination | approx. 110 ml  |
| Oil volume for dropping point determination            | approx. 100 ml  |
| Temperature sensor                                     | Pt-100  |
| Temperature range                                      | Room temperature up to 280° C   |
| Temperature program                                    | Rise of 0.5, 1, 2 or 3° C per minute;<br>30 to 40° C above the preset value |
| Regulating precision of set point potentiometer        | 0 ... 200° C: +/- 5° C<br>200 ... 280° C: +/- 7.5° C                        |
| Ramp precision   | +/- 10%   |
| Concordance thermometer - digital display              | 0...200° C: +/- 0.7° C<br>200...280° C: +/- 1° C                            |
| Cooling from 250° C to 60° C                           | approx. 13 minutes  |
| Heating from 60° C to 250° C                           | approx. 7 minutes   |
| Permissible ambient temperature                        | 0 ... 40° C   |
| Dimensions (width x height x depth)                    | 220 x 285 x 260 mm  |
| Weight   | 4.8 kg  |

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
**Instrument versions****BÜCHI 535**

For melting and boiling point determination

|       |          |        |
|-------|----------|--------|
| 220 V | 50/60 Hz | 22 470 |
| 240 V | 50/60 Hz | 22 471 |
| 117 V | 50/60 Hz | 22469  |

This instrument can be adapted at any time with the heating module for dropping point determination (see Accessories).

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
**2. Standard outfit (version for melting and boiling point determination)**

|        |  |   |
|--------|--|---|
| 1      | BÜCHI 535, completely assembled                          |   |
| 1      | Protective covering                                      | 25106                                     |
| 1      | Mains cable CH<br>or USA<br>or Schuko<br>or GB<br>or AUS | 10021<br>10023<br>10029<br>17833<br>17834 |
| 1      | Operating instruction, trilingual                        | 95340                                     |
| 1      | Holder for melting and boiling point tubes               | 20871                                     |
| 100    | Melting point tubes                                      | 17808                                     |
| 10     | Boiling point capillaries                                | 19697                                     |
| 10     | Boiling capillaries                                      | 19698                                     |
| 250 ml | Silicone oil   | 01757                                     |
| 2      | Spare fuses for mains voltage 220/240 V; 110 V           | 01687                                     |
| 1      | Antiglare device   | 22814                                     |
| 1      | Funnel   | 22899                                     |
| 1      | Light background screen                                  | 22796                                     |

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


### 3.1. Description of the controls

- 1 Main ON/OFF switch: The instrument is heated to the value preset by means of the dial (7), providing that the fan (2) is switched off.
- 2 When the fan is switched on the heater is simultaneously switched off.
- 3 Reset: The heating ramps (4) are switched off and the temperature drops to the preset value (7). The temperatures stored by means of the STO (5) button are retained.
- 4 Temperature program: This permits a temperature rise of 0.5, 1, 2 or 3°C per minute, up to 30-40°C above the preset value.
- 5 A maximum of three temperatures can be stored using this button.
- 6 RCL: By pressing the RECALL button, the temperatures stored using the STO button can be individually called up on the digital display.
- 7 Set-point potentiometer: Rapid heating to a preset value.
- 8 Heating indicator lamp:  
Lamp flashing = set value reached  
Lamp continuously lit = heating switched on, set value not yet reached  
Lamp not lit = heating switched off
- 10 Digital display: temperatures up to 325°C (280°C plus ramp of about 30°C) can be digitally displayed. If this temperature is exceeded, the display starts to flash.

#### 10a Indication of store - No 1-3

Erasing the stored values:

- a) After using the heating ramps (4): by pressing the O-button and then the STO button, the values stored in the memories 2 and 3 are erased. The value stored in the memory 1 is retained and overwritten with the instantaneous value.
- b) Without using the heating ramps beforehand: Press one of the heating ramp buttons (4), then the O-button and finally the STO button.

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#### 4. Installation of the B-530

##### 4.1. Instructions and special precautions

This instrument has been manufactured in accordance with the Safety Control Standards for Electronic Measuring Equipment, DIN 57411, Section 1, and VDE 0411, Section 1, and has left our factory in safe working order. To ensure that it remains in this condition, the user must observe the instructions and special precautions indicated below.

Before switching on the instrument, it must be ensured that the operating voltage for which the instrument has been set corresponds to the mains voltage.

The mains plug must only be plugged into a socket with a satisfactory earth connection. An extension cable without an earth lead must not be used.

##### Warning

Any break in the earth lead either inside or outside the instrument or a loose earth connection can lead to unsafe operation of the instrument. The earth lead must not be deliberately disconnected.

If it is necessary to remove the covers from the instrument in order to adjust, service, repair or replace parts, the instrument must first be disconnected from all power sources to avoid the danger of contact with live components. If it is required to service or repair an uncovered instrument which is still connected to a power source, then this should only be performed by a qualified engineer who is aware of the dangers.

It must be ensured that only fuses of the stated type and of the correct amperage are used as replacements. Repaired fuses should not be used, nor should short-circuiting of the fuse-holder be attempted.

If the instrument is no longer in safe working order, then the instrument must be taken out of service.

The instrument is no longer in safe working order if:

- it is visibly damaged;
- it no longer functions;
- it has been stored for a long period of time in adverse conditions;
- it has received heavy knocks during transport.

The instrument should be operated at temperatures above 200° only for short periods of time. This will keep the oil in good condition.

If the instrument is operated for any length of time at temperature above 200° C, the temperature safety switch may switch off the heating system and the motor.

Care should be taken when changing specimens, as there is the danger of coming into contact with hot surfaces.

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#### 4.2. Setting up the apparatus

Set up the apparatus in such a way that no light is able to fall directly onto the magnifying glass. This will ensure that there is no troublesome reflection from the magnifying glass.

Antiglare device: Loosen the screw on the antiglare device and place the screen over the magnifying glass (as shown in Figure). Tighten the screw again.

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#### 4.3. Changing the background screen

A black background screen is fitted as a standard feature. The white background screen is more suitable for dark substances.

Procedure:

- Remove the heating module as described in Section 11.2
- Release the spring grips from the glass vessel
- Release the wire arm 
- Pull out the heating vessel from underneath
- Pull out the dark screen from the top and insert the light screen 

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
#### **4.4. Filling the apparatus with silicone oil**

Insert the plastic funnel into the opening provided for the specimen holder.  
Twist the expansion vessel downwards into position.  
Pour silicone oil into the funnel until oil can be seen flowing into the expansion vessel.

#### **4.5. Mains connection**

Check that the mains voltage is the same as that indicated on the identification plate of the BÜCHI 535. Connect the B-535 to the means of the power supply cable (19)

### **5. Inserting the specimen tubes**

Insert the specimen tube inside the prism-shaped specimen holder so that it lies flush with the top edge of the metal part.

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#### **5.1. Inserting the specimen holder**

Insert the specimen holder vertically into the bracket and press it down fully into position.

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## 6.2. Determination of a substance with a known melting point

- Grind specimens which have large crystals and are non-homogeneous, in a mortar, and fill a specimen tube to a height of about 2 mm.
- Insert the specimen tube in the specimen holder.
- Completely empty the expansion vessel by briefly twisting it upwards.
- Set the set-point potentiometer (7) to a value about 10° C below the expected melting point and set the main switch to the ON position.
- Reset the temperature program by pressing the O- button (3).
- As soon as the oil temperature reaches the value selected on the set-point potentiometer, press the button for the required temperature rise; by pressing the RESET- button, it is possible to halt the temperature rise within a given temperature program; the temperature then drops back to the set value.
- Observe melting of the specimen and store the temperatures immediately at the beginning and end of the melting process by pressing the STO button (5).
- Then call up the measured temperatures by pressing the RCL button (6) and make a note of them.

### NB:

Using the temperature program, it is possible to achieve a temperaturerise of about 30° C. Switching on the cooling fan automatically switches off the heater.

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

### 6.3. Determining the melting point of a mixture

- Homogeneously mix in the mortar equal amounts of the two substances A and B whose melting point is to be determined.
- Fill three specimen tubes, each to a height of 2 mm, with the mixture, substance A and substance B. Insert the tubes into the specimen holder.
- Completely empty the expansion vessel by briefly twisting it upwards.
- Set the set-point potentiometer (7) to a value about 15° C below the expected melting point and set the main switch (1) to the ON position.
- Reset the temperature program by pressing the O- button (3).
- As soon as the oil temperature reaches the value selected on the set-point potentiometer, press the button for the required temperature rise and observe the specimen; by pressing the RESET- button, it is possible to halt the temperature rise within a given temperature program; the temperature then drops to the set value.
- Observe melting of all three specimens and store the temperatures immediately at the beginning and end of the melting process by pressing the STO button (5).
- Then call up the measured temperatures by pressing the RCL button (6) and make a note of them.

If the measured melting point of the mixture of A and B is significantly lower (at least 10° C lower) than the melting point of the lower melting, pure component, then A and B are two different substances. If the mixture has a melting point which is similar to or the same as that of the two components A and B, then the two substances are identical.

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
## 7. Performing a boiling point determination

- Fill the specimen tube (3.2 mm diameter) with the specimen fluid, to a height of 5-10 mm.
- Immerse the boiling capillary in the specimen fluid, thick end first.
- Insert the specimen tube in the specimen holder.
- Completely empty the expansion vessel by briefly twisting it upwards.
- Set the main switch (1) to the ON position and reset the temperature program by pressing the O-button(3).
- Set the oil temperature on the set-pointpotentiometer (7) to a temperature about 10° C below the expected boiling point. If the approximate boiling point is not known, a rapid preliminary determination is performed, as for melting point determination (Section 6.1).
- When the set value is reached, press the program button (4) for the required temperature rise and observe the specimen.
- The boiling point is reached when there is a fast continuous stream of bubble (like a pearl necklace).
- Store the corresponding temperature by pressing the STO button (5).
- Then call up the temperature measured by pressing the RCL button (6) and make a note of it.

### NB:

Using the temperature program, it is possible to achieve a temperature rise of about 30-40° C. Switching on the cooling fan automatically switches off the heating system.

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
## **8. Performing a dropping point determination**

(using the Ubbelohde method, DIN 51801, Section 2)

Insert the heating module for dropping point determination, as described in Section 11.2.

NB: When the heating module for dropping point determination is inserted, the digital display no longer functions.

### **Preparing the specimen:**

Place the nipple on a smooth glass or metal surface and fill it with the specimen previously heated to the liquid state.

NB: In the case of sticky specimens, coat the plate beforehand with a release agent (for example a 1:1 mixture of glycerine and dextrin).

Allow the filled nipple to cool and place it for about 30 minutes in a refrigerator or in ice-water.

Using a heated knife, trim off any excess specimen so that the surface of the specimen is smooth.

At a suitable temperature, push the nipple onto the holder located on the thermometer.

Insert the entire thermometer assembly into the glass vessel.

Set the main switch (1) to the ON position and reset the temperature program by pressing the O- button (3).

Set the set-point potentiometer (7) to a temperature about 10° C below the expected dropping point.

As soon as the set temperature has been reached, press the 1° C/min. temperature program button.

The dropping point is the temperature shown on the thermometer when the first drop falls from the nipple opening.

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## **9. Maintenance and cleaning**

### **9.1. General**

The housing of the Büchi 535 used for determining melting points has a painted finish.

Cleaning: Use only a dry cloth or alcohol.

### **9.2. Changing the oil**

Silicone oil can with time polymerise, especially if the instrument is operated at high temperatures. This increases the viscosity of the oil and hence the oil must be periodically changed.

#### **Procedure:**

- Pull out the mains plug.
- Remove the heating module as described in Section 11.2.
- Remove the expansion vessel by releasing the clips from the ground-glass joint.
- Remove the oil via the overflow opening.
- Release the three spring grips on the glass vessel and swing back the wire arm.
- Hold the glass vessel and now pull out the motor from the top and place it in a suitable beaker.

#### **Cleaning of the glass vessel:**

Rinse out the heating vessel using a mixture of petroleum ether and methylene chloride.

### **9.3. Cleaning the magnifying glass**

The magnifying glass of an acrylic lens and is not resistant to solvents.

Use only ethanol for cleaning.

**Do not use acetone, chlorinated solvents or similar agents.**

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M1 Stirrer motor  
M2 Fan  
E1 Heating  
E2 + E3 Illumination  
R45 Regulating sensor  
R46 Measuring sensor  
S20 Overheating safety switch  
X4 Plug illumination and fan  
X5 Plug for digital section  
X7 Plug for display PCB

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**Colour code for cables**

|              |
|--------------|
| yellow/green |
| brown        |
| blue         |
| black        |
| white        |
| red          |
| green        |
| grey         |
| violet       |
| yellow       |
| pink         |

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
|       |                      | 117/220/240 V | Set of ten |
|-------|----------------------|---------------|------------|
| F1/F7 | Main fuse            | 2.5 AT        | 22560      |
| F2    | Heater               | 2.5 AFF       | 22570      |
| F3/F4 | Analog Section       | 0.1 AT        | 22555      |
| F5    | Digital Section      | 0.5 AT        |            |
| F6    | Fan and illumination | 0.5 AT        | 22558      |

|       |                         | 117/220/240V | 10 pièces |
|-------|-------------------------|--------------|-----------|
| F1/F7 | Fusible principal       | 2.5 AT       | 22560     |
| F2    | Chauffage               | 2.5 AFF      | 22570     |
| F3/F4 | Partie analogique       | 0.1 AT       | 22555     |
| F5    | Partie digitale         | 0.5 AT       |           |
| F6    | Ventilateur + Eclairage | 0.5 AT       | 22558     |

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

## **11. Replacing parts**

### **11.1. General precautions**

Before replacing any parts on the instrument, always:

- Set the mains switch to the OFF position
- Pull out the mains plug

### **11.2. Changing the heating module**

Pull out the temperature sensor plug (17). Using a screwdriver, loosen the screw on the cable (16) for the motor and heater and pull out the plug. Release the quick fasteners (o) by applying light pressure and at the same time performing a half screwdriver turn. Pull out the heating module from the top.

### **11.3. Replacing the magnifying lens**

- Loosen the screws on the front plate on both sides.
- Swing out the frame with the magnifying lens and the digital display.
- Loosen the screws on the magnifying glass mounting.
- Remove the magnifying glass with its holder.
- Change the magnifying glass.

### **11.4. Replacing the cooling fan**

- Remove the heating module by releasing the Camlock fasteners.
- Pull out the flat push-on sleeves on the fan.
- Loosen the screws (o).
- Mount the new fan.

**Note:** If the flat plug is not mounted properly, the circuit board may be damaged.

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### **11.5. Replacing the heating coil**

Remove the heating module as described in Section 11.2.

Remove the spring grips and the wire arm from the heating vessel.

Pull out the glass vessel from underneath.

Press out the heating coil from the mounting with the help of a small screwdriver.

Press new heating carefully into holder.

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### **11.6. Replacing the illumination lamps**

Remove the heating module, as described in Section 11.2.

The lamp and reflector are of the push-in type. Pull the lamp with reflector out and put a new one in its place. It is also possible to change the lamp by opening the front frame (see Section 11.3).

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### **11.7. Replacing the Digital display PCB GS 3**

- Loosen the screws on the front plate with the magnifying glass and digital display
- Swing out the front plate
- Loosen the screws on both sides of the housing (o) and lift off the top section of the housing
- Unscrew the digital PCB (o)
- Disconnect the connector from the main PCB; to do this, press both retainers outwards
- Unscrew and remove the metal plate (o)
- Pull out the connector through the opening and replace PCB

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## 12. Troubleshooting

| Problem   | Possible cause   | Solution  |
|---|--|---|
| Instrument not working at all                                     | 1.No current<br>2.Fuse F1/F7 defective   | 1. Check mains voltage<br>2. Change fuse  |
| Stirrer and illumination working,<br>heating not working          | 1.Fuse F2, F3 or F4 defective<br>2.Heating coil defective<br>3.Electronics defective<br>4.Pt-100 defective | 1. Change fuse F2, F3 or F4<br>2. Replace heating coil<br>(see Section 11.5.)<br>3. Call Service Engineer<br>4. Call Service Engineer |
| Illumination not working  | 1.Lamps defective<br>2.Fuse F6 defective<br>3.Plug X4 not connected to PCB                                 | 1. Change lamps (see Section 11.6)<br>2. Change fuse<br>3. Connect plug to PCB  |
| Stirrer not working,<br>heating working                           | Stirrer motor defective  | Call Service Engineer   |
| Fan not working   | 1.Fuse F6 defective<br>2.Fan defective<br>3.Plug X4 not connected to PCB                                   | 1. Change fuse<br>2. Replace fan (see Section 11.4)<br>3. Connect plug to PCB   |
| Set value incorrect   |  | Call Service Engineer   |
| Temperature heating rate incorrect                                |  | Call Service Engineer   |
| Stirrer and heating not working                                   | 1.Cable 16 not connected<br>2.Overheating safety switch<br>has come into operation                         | 1. Connect cable<br>2. Allow instrument to cool   |
| Digital display not working                                       | 1. Flat cable not connected<br>2. Digital display PCB defective<br>3. Fuse F5 defective                    | 1. Connect cable X5 Or X7 (see Section 10.1.)<br>2. Replace PCB (see Section 11.7.)<br>3. Change fuse F5                              |
| Flashing error code (000/399) appears<br>on digital display       | 1. Cable 17 not connected<br>2. Pt-100 defective   | 1. Connect cable<br>2. Call Service Engineer  |
| Incomplete information shown<br>on digital display                | LED defective  | Replace PCB (see Section 11.7.)   |
| Different values obtained from<br>digital display and thermometer | 1. Oil to viscous<br>2. Thermometer imprecise<br>3. Electronics disadjusted                                | 1. Change oil (See Section 9.2.)<br>2. Use calibrated thermometer<br>3. Call Service Engineer   |
| Temperatur not stable<br>(Indication raises and drops)            | Circulation of the oil is<br>not sufficient  | 1. Adjust stirrer in the heating vessel<br>2. Change oil (viscosity too high)   |

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